

CASO 12:



Operadores: **Raquel del Valle Fernández**
Pilar Jiménez Quevedo

Moderadores:

- **Fernando Sarnago Cebada**
- **Bruno García del Blanco**

Panelistas:

- **Pablo Avanzas Fernandez**
- **Roberto Blanco**
- **María Cruz Ferrer García**
- **Cristobal Urbano Carrillo**
- **Rafael Rodríguez de la Borbolla**

Glidesheath Slender® Hydrophilic Coated Introducer Sheath

- 1 **Diámetro externo reducido en 1 Fr**
- 2 **Mantiene el diámetro del lumen**
- 3 **Reduce la necesidad de cambiar a un introductor de talla superior**

Introductor
convencional
6Fr

**Lumen
Equivalente**

Glidesheath
Slender®
6Fr

**Diámetro externo
Equivalente**

Introductor
convencional
5Fr


SAFARI²™

Pre-Shaped Guidewire


Enhanced wire predictability
with superior shape retention

Streamlined device delivery
through optimized rail support

Widest guidewire
choice with three curve sizes



- Outer Diameter: 0.035" (0.889 mm)
- Overall Length: 275 cm
- Core Material: Stainless Steel
- Coil Material: Stainless Steel
- Coating: LUBRIGREEN™ PTFE
- Unique Product Dimensions by Curve Size





**Boston
Scientific**

Heart- Navigator

Fusing live X-ray and live

PHILIPS Healthcare





-  High pressure and high resistance balloon
-  Specially designed for pre-dilation of the native valve prior to prosthesis implantation
-  Advanced folding with flat folds for easy removal after dilation
-  Fast inflation and deflation times
-  Soft catheter tip and short sheath

Evolut™ FX+ TAVI System

5, 6 y 7 NOVIEMBRE
HOTEL RIU PLAZA DE ESPAÑA

Medtronic

Outflow

Designed with the
same structural strength
at the outflow as the
Evolut™ FX valve.¹

Built for durability

Inflow

Designed with the **same
radial force** at the inflow
as the Evolut™ FX valve.¹



Designed for access

Waist

Designed with three
windows that are
4x larger to enable
lifetime management
solutions, such as
coronary access.²

Evolut™ FX+ TAVI System



Evolut FX Loading System

Model number:
L-EVOLUTFX-2329


23, 26, and 29 mm TAV



Evolut FX Loading System

Model number:
L-EVOLUTFX-34

34 mm TAV



Medtronic

Evolut™ FX+ TAVI System

Medtronic



Evolut FX+ 23 mm Valve

Model number:
EVFXPLUS-23

Size: 23 mm

Aortic annulus
diameter:
17⁺/18-20 mm



Evolut FX+ 26 mm Valve

Model number:
EVFXPLUS-26

Size: 26 mm

Aortic annulus
diameter:
20-23 mm



Evolut FX+ 29 mm Valve

Model number:
EVFXPLUS-29

Size: 29 mm

Aortic annulus
diameter:
23-26 mm



Evolut FX+ 34 mm Valve

Model number:
EVFXPLUS-34

Size: 34 mm

Aortic annulus
diameter:
26-30 mm

BUILT UPON 29+ YEARS OF THE PERCLOSE™ LEGACY^{1*}

Perclose™ ProStyle™ Suture-Mediated Closure and Repair System



HIGHER TENSILE STRENGTH NEEDLES

15% increased suture deployment success
in the same complex, challenging anatomy¹

IMPROVED EASE OF USE

Improved usability and more intuitive
deployment experience¹

SHORT PATIENT RECOVERY

As short as 2-hour time to ambulation for small
and
large-bore procedures²

15
MILLION
REPAIRS^{1**}

